

SAMO ZA STRUČNU UPORABU



# AMINAL Cart

Kapsula s natrijevim bikarbonatom -  $\text{NaHCO}_3$

HEMODIJALIZA

# AMINAL Cart

Kapsula s natrijevim bikarbonatom - NaHCO<sub>3</sub>

## TEHNIČKI PODACI

Opis	AMINAL Cart kapsula s natrijevim bikarbonatom za hemodijalizu									
	620 g	650 g	720 g	750 g	760 g	900 g	950 g	1000 g	1100 g	1150 g
Kataloški broj	CART.050	CART.060	CART.057	CART.052	CART.053	CART.058	CART.055	CART.054	CART.051	CART.059
Sastav	Natrijev bikarbonat u prahu NaHCO <sub>3</sub>									
Namjena	On-line priprema tekućeg koncentrata bikarbonata za hemodijalizu									
Pakiranje - kutija koja sadrži	10 kaps.	10 kaps.	10 kaps.	10 kaps.	10 kaps.	8 kaps.	8 kaps.	8 kaps.	8 kaps.	8 kaps.
Kompatibilan s	Gambro [AK10, AK200, AK200S, AK200S ULTRA, AK90, AK95, AK96, AK98, Phoenix, Artis Physio], Hospal [Integra], B.Braun [Dialog+, Dialog IQ], Bellco [Formula, Flexya (tvornička opcija)], Nikkiso [DBB-03, DBB-05, DBB-06, DBB-07, DBB-EXA], Nipro [Surdial 55plus (tvornička opcija), Surdial X, NCU-18, DIAMAX (tvornička opcija)], Dialife [Dianova, Dianova X2]					Gambro [AK10, AK200, AK200S, AK200S ULTRA, AK96, AK98, Phoenix, Artis Physio], B.Braun [Dialog+, Dialog IQ], Bellco [Formula, Flexya (tvornička opcija)], Nikkiso [DBB-03, DBB-05, DBB-06, DBB-07, DBB-EXA], Nipro [Surdial 55plus (tvornička opcija), DIAMAX (tvornička opcija)], Dialife [Dianova, Dianova X2]				
Rok uporabe	36 mjeseci									
Uvjeti čuvanja	5°C - 30°C									

Proizvod je u skladu sa zahtjevima Europske farmakopeje (Ph. Eur) i Američke farmakopeje (USP).



ALKALOID  
SKOPJE

*Health above all*

[www.alkaloid.com.mk](http://www.alkaloid.com.mk)

CE 2803

Verzija 01 iz 10.2023.

PHARMACEUTICAL CHEMICAL COSMETIC INDUSTRY  
Blvd. Aleksandar Makedonski 12, Skopje Republic of North Macedonia  
Tel.: +389 2 31 04 041; Fax.: +389 2 31 04 041