

FOR PROFESSIONAL USE ONLY



## AMINAL Cleaner A

**CARTRIDGE  
CONTAINING POWDER  
OF ANHYDROUS  
SODIUM CARBONATE  
– Na<sub>2</sub>CO<sub>3</sub>**

FOR PROFESSIONAL USE ONLY



## AMINAL Cleaner C

**CARTRIDGE  
CONTAINING POWDER  
OF ANHYDROUS  
CITRIC ACID  
– C<sub>6</sub>H<sub>8</sub>O<sub>7</sub>**

# HAEMODIALYSIS

# AMINAL Cleaner A

**CARTRIDGE CONTAINING  
POWDER OF ANHYDROUS  
SODIUM CARBONATE**

**– Na<sub>2</sub>CO<sub>3</sub>**

# AMINAL Cleaner C

**CARTRIDGE CONTAINING  
POWDER OF ANHYDROUS  
CITRIC ACID**

**– C<sub>6</sub>H<sub>8</sub>O<sub>7</sub>**

## TECHNICAL DATA - AMINAL Cleaners

Description	AMINAL Cleaner A Cartridge containing powder of anhydrous sodium carbonate – Na <sub>2</sub> CO <sub>3</sub>	AMINAL Cleaner C Cartridge containing powder of anhydrous citric acid – C <sub>6</sub> H <sub>8</sub> O <sub>7</sub>
<b>Composition</b>	Sodium Carbonate (Na <sub>2</sub> CO <sub>3</sub> )	Citric Acid (C <sub>6</sub> H <sub>8</sub> O <sub>7</sub> )
<b>Model no.</b>	CART.CL51	CART.CL52
<b>Intended use</b>	On-line production of sodium carbonate solution intended for the disinfection and cleaning of dialysis machines, as well as for the removal of organic deposits from the fluid path of the dialysis machines.	On-line production of citric acid solution intended for the disinfection and cleaning of dialysis machines, as well as for the removal of calcium and magnesium salts from the fluid path of the dialysis machines.
<b>Packaging</b>	13 g	32 g
<b>Compatible with</b>	Gambro dialysis machines	
<b>Shelf life</b>	2 years	
<b>Storage conditions</b>	5°C - 30°C	

The product is in compliance with the requirements of European Pharmacopoeia and USP.



**ALKALOID  
SKOPJE**

*Health above all*

[www.alkaloid.com.mk](http://www.alkaloid.com.mk)

**CE 2803**

Version 03 from 12.2023.

PHARMACEUTICAL CHEMICAL COSMETIC INDUSTRY  
Blvd. Aleksandar Makedonski 12, Skopje Republic of North Macedonia  
Tel.: +389 2 31 04 041; Fax.: +389 2 31 04 041